

Clinical studies confirm the difference


More than 69 studies with over 3,000 test subjects scientifically confirm what users and technicians have experienced with more than 75,000 C-Leg fittings around the world.

Further information:
www.ottobock.com/clinicalstudies/lower-limb-prosthetics



References

1. Wismer, N., Mileusnic, M., Sreckovic, I. & Hahn, A. (2016). First results on next generation C-Leg. Poster presentation at OTWorld, Leipzig, Germany.
2. Kahle, J. T., Highsmith, M. J. & Hubbard, S. L. (2008). Comparison of nonmicro-processor knee mechanism versus C-Leg on Prosthesis Evaluation Questionnaire, stumbles, falls, walking tests, stair descent, and knee preference. Journal of Rehabilitation Research and Development, 45(1), 1-14.
3. Drerup, B., Wetz, H. H., Bitterle, K. & Schmidt, R. (2008). Langzeitergebnisse mit dem C-Leg – Ergebnisse einer retrospektiven Studie [Long-Term Results with the C-Leg – Results of a Retrospective Study] Orthopädie-Technik, 3, 169-174.



4. Highsmith, M. J., Kahle, J. T., Bongiorno, D. R., Sutton, B. S., Groer, S., Kaufman, K. R. (2010). Safety, energy efficiency, and cost efficacy of the C-Leg for transfemoral amputees: A review of the literature. Prosthetics and Orthotics International 2010; 34(4): 362-377.

ottobock.

C-Leg 4 – Reclaim your determination

Know the difference

TECHNICAL UPDATE 2018



Information for practitioners

Benefits

- ▶ Reliable stance and swing phase control in real time on various surfaces
- ▶ Choice between Intuitive and Manual Stance
- ▶ Safe backward walking
- ▶ Active Stumble Recovery at all times
- ▶ Cockpit app for Android and iOS devices*
- ▶ C-Soft Plus adjustment software with video tutorials supports the fitting process (version 1.6 and up)
- ▶ Large selection of compatible feet
- ▶ Weatherproof with IP67
- ▶ Two MyModes



New C-Leg 4: What you need to know

- Product reference numbers (3C98-3*, 3C88-3*) remain the same
- No recertification necessary
- Update of C-Soft Plus to version 1.6 (4X440=V1.6) required
- New Cockpit app available in the Google Play Store (Android) and Apple App Store (iOS)*
- Remote control (4X350) cannot be used for iOS-compatible C-Leg

For more information, see the instructions for use or the new C-Soft Plus version 1.6.

*Apple iOS compatibility is available with the new C-Leg 4 only. iOS compatibility is not available with previous versions of the C-Leg 4.

Safety first:
More robust design for better protection against external influences.

New battery covers:
Even better protection for the valuable interior.

Cockpit app:
Now available for Android and Apple iOS devices.

Video tutorials in the new C-Soft Plus:
Clear explanation of the fitting process.

Even larger selection of prosthetic feet:
19 predefined feet (a total of 27 versions) for a faster, individualized fitting.



Shield insert
Art. no. 4P863=*

C-Leg Guard | Art. no. 4P862
Guard against external influences.



Shield insert
Art. no. 4P863=*

C-Leg Protective Cover | Art. no. 4X860=*
Individually adaptable protector with foot cuff for C-Leg and tube adapter.

